

Shugga for Loop

Shugga for Loop is a simple yet highly customizable iPhone app that announces your latest blood glucose (BG) data found in your HealthKit. Primarily built as a companion for Loop.app, it also serves those who have their live BG values written to HealthKit by other means, and even caters to those who manually enter their BG values into Health app, though in a limited capacity.

Additionally, it provides an optional display of the carbohydrates you've ingested in the past 24 hours, also sourced from Health data.

The app's design and customization cater to a wide range of users, offering enhanced accessibility for visually impaired individuals and features for active users. The main page is deliberately kept simple for quick visual cues.

Shugga for Loop is exclusively compatible with iPhones due to its dependence on Apple Health, which is not available on other devices.

Shugga for Loop is an independent app and is not affiliated with Loop.app, LoopKit, or Dexcom.

We take your privacy seriously and want you to know that Shugga for Loop does not use your data for any other purpose. Your data is never stored, transmitted, or shared outside the confines of this app. Your privacy is important to us and we are committed to keeping your data to you only.

You will need to grant the app an ongoing "Read Access" to your BG data in HealthKit for this app to work at all. Optionally, you can give the same access to your carbohydrate data if you are interested in being able to see your carbohydrate history in the app.

You will be prompted only once for the permissions by the Shugga for Loop, per Apple's design at the very first time you use the app. If you don't grant access(es) then, you can do so later in the iPhone's Settings: Health: Data Access and Devices: Shugga for Loop.

If you want Shugga for Loop to work while the app is in the background, make sure you enable Background App Refresh in the iPhone's Settings: Shugga for Loop: